

## **Device Use Case : Cranial Electrotherapy Device**

### **Turning Device On (UC1)**

Primary Actor:

- Device User

Precondition:

- User has Device

Main success scenario:

1. User turns on the device by pressing “On” Button.
2. The device loads up the display elements and displays information to user.

Postcondition:

- Device has been turned on and waits for user input to start treatment.

Extensions:

- 2a. Device does not turn on.
  - 2a1. Charge device to ensure the battery level is above 2%.
  - 2b. Device detects current overload.
    - 2b1. Device permanently disables device.
    - 2b2. Contact manufacturer to repair device.

Related Information:

- Device will show a warning when the battery is less then 5% stating that “ The battery is low”

### **Turning Device Off (UC2)**

Primary Actor:

- Device User

Precondition:

- Device is turned on

Main success scenario:

1. User turns off the device by pressing “Off” Button.
2. The device shutdown and does not display information anymore.

Postcondition:

- Device has been turned off and doesn’t consume power.

Related Information:

- Device will turn of after 30 minutes of being not in use.

### **Changing Treatment Time On Device (UC3)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device and is in the main menu.

Main success scenario:

1. User selects the wanted duration of therapy.
2. Device updates the view to new duration.
3. Device sets the settings to the new duration.

Postcondition:

- New settings now set to new duration of therapy.

Related Information:

- Treatment of time can only be set once the treatment is finish or cancelled.

### **Changing Frequency On Device (UC4)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device and is in the main menu.

Main success scenario:

1. User selects the wanted frequency.
2. Device updates the view to new duration.
3. Device sets the settings to the new duration.

Postcondition:

- New setting is now set to the new frequency.

Related Information:

- Frequency can be change anytime during the treatment.

### **Changing Waveforms On Device (UC5)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device and is in the main menu.

Main success scenario:

1. User selects the wanted waveforms.
2. Device updates the view to the new waveform.
3. Device sets the settings to the new waveform.

Postcondition:

- New setting is now set to new waveform.

Related Information:

- waveform can be change anytime during the treatment.

### **Decreasing Intensity on Device (UC6)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device and is in the main menu.

Main success scenario:

1. User presses "Decrease Intensity".
2. Device will decrease the intensity of the device by 1.
3. Device will use the new setting in the treatment.

Postcondition:

- The device has stored the new intensity in the settings and will run using it when treatment begins/while running.

Extensions:

2a. Intensity does not decrease

2a1. The intensity of the device is already the lowest, zero.

Related Information:

- The "Decrease Intensity" button is a cycle, pressing it once will decrease it by one until it is 0 where it will no longer go down and the button disables.
- Each intensity is 100 microampere on top of the previous one e.g., Intensity 3 is 300 microampere
- Intensity can be change at any time during the treatment.

### **Increasing Intensity On Device (UC7)**

**Primary Actor:**

- Device User

**Precondition:**

- User has turn on the device and is in the main menu.

**Main success scenario:**

1. User presses “Increase Intensity”.
2. Device will increase the intensity of the device by 1.
3. Device will use the new setting in the treatment.

**Postcondition:**

- The device has stored the new intensity in the settings and will run using it when treatment begins/while running.

**Extensions:**

2a. Intensity does not increase

2a1. The intensity of the device is already the highest , five.

**Related Information:**

- The “Increase Intensity” button is a cycle, pressing it once will increase it by one until it is 5 where it will no longer go up and the button disables.
- Each intensity is 100 microampere on top of the previous one e.g., Intensity 3 is 300 microampere
- Intensity can be change at any time during the treatment.

### **Viewing Treatment History On Device (UC8)**

**Primary Actor:**

- Device User

**Precondition:**

- User has turn on the device and is in the main menu.

**Main success scenario:**

1. User presses “View History”.
2. Device displays a list of all the history that ran on the device.

**Postcondition:**

- The device displays the stored history of all the treatments that was ran.

**Related Information:**

- Can be viewed after or during the treatment.

### **Returning to Main Menu from Treatment History Page On Device (UC9)**

**Primary Actor:**

- Device User

**Precondition:**

- User is in the “View History” display.

**Main success scenario:**

1. User presses the “Back” button.
2. Device returns user to the main window.

**Postcondition:**

- Returns to main window.

**Related Information:**

- Back button is only used for returning from “view history” window.

### **Locking Settings Buttons On Device (UC10)**

**Primary Actor:**

- Device user

Precondition:

- User is running a treatment.

Main success scenario:

1. User presses the “lock” button.
2. Device locks the device from any settings change.

Postcondition:

- Settings buttons become unresponsive.

Related Information:

- Lock is only the buttons that effect changes in the settings of the treatment. - Pressing it again will unlock the settings.

### **Record Treatment On Device (UC11)**

Primary Actor:

- Device user

Precondition:

- User is running a treatment.

Main success scenario:

1. User presses the “Record” button.
2. Device stores all the settings of the treatment in the database.

Postcondition:

- Settings are stored for later viewing.

Related Information:

- Only records the settings that were set before pressing record, further changes after record is pressed does not get recorded in the history.

### **Applying Nodes to User’s Skin (UC12)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device and is in the main menu.

Main success scenario:

1. Device user attaches electrodes to skin with proper contact.

Postcondition:

- Device detects that proper skin contact is achieved.

Extensions:

- 1a. Attach electrodes do not make proper contact.
  - 1a1. Remove electrodes from skin
  - 1a2. Clean surface of skin from any debris.
  - 1a3. Apply pre-package oil to target area to ensure proper contact to skin.

Related Information:

- Treatment will only start when proper contact is made and intensity are set.

### **Removing Nodes From User’s Skin (UC13)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device and is in the main menu.

Main success scenario:

1. Device user detaches electrodes from skin.

Postcondition:

- Device detects contact have been removed.

Extensions:

- 1a. Electrodes are detached from skin for longer than 5 seconds.
  - 1a1. Device detects lost of contact to skin
  - 1a2. Device terminates treatment and doesn't send current through.

Related Information:

- If the electrodes are removed for less than 5 seconds, current will still run through it.

### **Device Initiates Treatment (UC14)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device, attach the electrodes properly and set intensity to higher than 0.

Main success scenario:

1. Device starts treatment.

Postcondition:

- Device runs current through electrodes with the settings that were set.

### **Device Terminates Treatment (UC15)**

Primary Actor:

- Device User

Precondition:

- User has turn on the device, set intensity to higher than 0 and detaches the electrodes while the treatment is running.

Main success scenario:

1. Device stops treatments.

Postcondition:

- Device stops sending current through the electrodes.

Relevant Information:

- Reattaching the electrodes will continue the treatment.

### **Device Update The Display (UC16)**

Primary Actor:

- Device User

Precondition:

- User changes settings or activates a button on main menu while turn on with sufficient battery.

Main success scenario:

1. Device will update the appropriate sections with the newly selection option or settings

Postcondition:

- Display shows new information or request.

